

- Patented item
- A uniquely designed interdental brush with the ability to be used at different angles to make cleaning easier.
- With a replaceable brush head making it more economical and environmentally friendly.
- Complete with 5 brush heads and a smart travel case. Ideal for travelling.





IB306

- The push up brush is used for cleaning between the teeth, to stimulate gums and effectively remove plaques.
- The push down dental pick is used to pick away annoying food particles between the teeth.





To remove unwanted objects from between the teeth, slowly insert the interdental brush in



Gemtly brush between all your teeth several times to remove the food and plaque and massage gums.



Rinse it off in the sink and push back in place on cover the cap.



www.caradent.com TEL : 886-7-555 9928 FAX : 886-7-555 9926